

# ORTHOFLEX tendon<sup>®</sup>

- ORTHOFLEX tendon<sup>®</sup> is viscoelastic gel for peritendinous or intrasheath injections.
- ORTHOFLEX tendon<sup>®</sup> is injectable prefilled syringe 40mg/2ml natural solution of highly purified sodium hyaluronate obtained by biofermentation devoid of animal protein + 10 mg mannitol, a free radical scavenger which help to stabilize the chains of sodium hyaluronate.

orto brand  
enjoy your life



[www.ortobrand.com](http://www.ortobrand.com)

## Indications:

- For the treatment of pain and restricted mobility in tendon disorders.
- Ankle and Foot
  - achilles tendinopathy
  - posterior tibial tendinopathy
  - peroneal tendinopathy
- Shoulder- rotator cuff tendinopathy
- Elbow- lateral and medial epicondylalgia
- Knee – patellar tendinopathy

ORTHOFLEX tendon<sup>®</sup> should be administered **once a week** for total of **2 injections**.

## Effects:

- Acts a **lubricant** when applied into the tendon sheath or peritendinously.
- Enhance tendon **gliding effect** and **reduce adhesion**.
- **Prevents** the free passage of **inflammatory cell** and molecules through the tendon sheath due to its macromolecular meshwork.
- Promotes tendon recovery and wound healing process.

## Precautions:

ORTHOFLEX tendon should be instilled accurately into tendon sheath or around the affected tendon if necessary under imaging control, avoid nerve lesions and injections into blood vessels.

## Characteristics:

- Treatment of tendon disorders
- **2 ml prefilled syringe for single use**
- Sodium hyaluronate 2% (40mg/2ml) + mannitol (10mg)
- Molecular weight- 3000kDa
- **Duration** on effects – **6 month**

## Quality:

ORTHOFLEX tendon<sup>®</sup> is produced according to Good Manufacturing Practice (GMP) and ISO13485 being certified as a system with high quality ingredients, CE marked medical device being in accordance with the European Directive 93/42/EE